

What's  
new?

## 26" Full HD LCD Monitor

Ackermann®



### FEATURES

- ▶ Calibrated to clinical colour
- ▶ Image pan, zoom, freeze, picture-in-picture
- ▶ Rapid signal detection, robust mode tables
- ▶ Built to last, rugged aluminum housing
- ▶ Artifact-free HD images
- ▶ Fanless - sterile field compatible

### ADVANTAGES IN COMPARISON TO 16-2033HD-26

- ▶ Higher intensity of black colour
- ▶ More powerful colour depth
- ▶ Considerably lower price

### TECHNICAL DATA

|                        |   |
|------------------------|---|
| <b>Article number</b>  | <b>16-2033HD-26NT</b>   |
| Screen                 | 26"   |
| Resolution             | 1920 x 1080 pixel (HD 1080)   |
| Aspect ratio           | 16:9  |
| Active area            | 576.0 X 324.0 mm  |
| GPIO port              | Swap, PIP/ PBP1/ PBP2 Select, Record indicator  |
| Output signal          | 1 x DVI-D, 1 x SD/ HD/ 3G-SDI (BNC)   |
| Control key            | Power, Menu, PIP, Up, Down, Plus, Minus, Input  |
| Dimensions             | 638 (W) x 389 (H) x 74.7(D) mm<br>25.12(W) x 15.32(H) x 2.94(D) inch  |
| Input signal           | 1 x DVI-D 1<br>1 x DVI-D 2 (optical fiber optional)<br>1 x VGA (D-sub)<br>1 x SD/HD/3G-SDI (BNC)<br>1 x C-Video (BNC)<br>2 x S-Video (Y/C) (BNC)<br>1 x Component (R/GBS, Y/PbPr) (5 x BNC) |
| Average brightness     | 500 (5 points avg.) cd/m <sup>2</sup>   |
| Response time (G-to-G) | 8 ms (avg.)   |
| Viewing angle (CR>10)  | R/L 178, U/D 178  |
| Weight                 | 7.6 kg, 16.8 lbs  |
| Contrast ratio         | 1500 : 1 typical  |
| Pixel pitch            | 0.3 x 0.3   |
| Power supply           | AC/DC Adaptor (AC 100~240V, DC 24V 6.25A)   |
| Optional               | Footstand 16-2033HD-26FUSS  |

### Ackermann Instrumente GmbH

Eisenbahnstraße 65-67  
78604 Rietheim-Weilheim  
Germany

### Contact details

Phone +49 (0)7461 96617 0  
Fax +49 (0)7461 96617 70  
E-Mail info@ackermanninstrumente.de