

SHAVER SYSTEM & BLADES

INVESTMENT OPTIONS | SHAVER SYSTEM

3 Investment Options - All Opportunities

Option 1

- ▶ Initial Investment, Owner of Unit

Option 2

- ▶ No Investment, Owner of Unit

Option 3

- ▶ Investment of Installments, Only Loan Unit

16-2050

Shaver System incl. Handpiece & Footswitch

Facts & Features

- ▶ Speed adjustable from 800 – 8000 RPM
- ▶ Ergonomic and light weight handpiece
- ▶ Silent and smooth motion
- ▶ Handpiece (with cable), autoclavable
- ▶ All functions can be either commanded on the control unit
- ▶ Shaver offers the possibility to connect 2 tools at the same time (for example handpiece and shaver blade)



Technical Data

Main voltage	230 VAC
Power consumption	50 VA
Frequency	50 - 60 Hz
Fuses	2 x 1,6 AT
Dimensions	90 x 360 x 268 mm
Operation mode	continous operation
Connections	2x handpiece, 1x footswitch
Weight	6 kg
Device type	B
Casing safety standard	IPX0
Safety class	Ila
Power cord	230 V, 10 A, 3-pole

Included

- ▶ Handpiece (16-2050-100)
- ▶ Footswitch (16-2050-200)
- ▶ Autoclavable silicone tubing
- ▶ Power cord
- ▶ Fuses
- ▶ Instruction for use

SHAVER SYSTEM & BLADES

BLADES

Full Radius Recetor

Plain sheath tip, plain inner tip



Reusable

15-1722-100ART

15-1722-101ART

15-1722-102ART

Disposable

15-1722-100DDART

15-1722-101DDART

15-1722-102DDART

Dimensions

5.5 mm

4.5 mm

3.5 mm

Colour Code



Meniscus Cutter

Serrated sheath tip, plain inner tip



Reusable

15-1722-200ART

15-1722-201ART

15-1722-202ART

Disposable

15-1722-200DDART

15-1722-201DDART

15-1722-202DDART

Dimensions

5.5 mm

4.5 mm

3.5 mm

Colour Code




SHAVER SYSTEM & BLADES

BLADES

Aggressive Full Radius Resector

Plain sheath tip, serrated inner tip




Reusable	Disposable	Dimensions	Colour Code
15-1722-300ART	15-1722-300DDART	5.5 mm	
15-1722-301ART	15-1722-301DDART	4.5 mm	
15-1722-302ART	15-1722-302DDART	3.5 mm	

Incisor

Serrated sheath tip, serrated inner tip



Reusable	Disposable	Dimensions	Colour Code
15-1722-400ART	15-1722-400DDART	5.5 mm	
15-1722-401ART	15-1722-401DDART	4.5 mm	
15-1722-402ART	15-1722-402DDART	3.5 mm	

SHAVER SYSTEM & BLADES

BLADES

Round Burr Drill



Reusable

15-1722-500ART

Disposable

15-1722-500DDART

Dimensions

5.5 mm

Colour Code



Tapered Burr Drill



Reusable

15-1722-600ART

Disposable

15-1722-600DDART

Dimensions

5.5 mm

Colour Code



Oval Burr Drill



Reusable

15-1722-700ART

Disposable

15-1722-700DDART

Dimensions

5.5 mm

Colour Code

