


XPRESS LOCK ©

DETACHABLE INSTRUMENT LINE

new!



WHY XPRESS LOCK © ?

- › innovative quick-locking mechanism
- › time-saving - fast cleaning & assembly
- › short, standard and bariatric lengths
- › no hidden space to accommodate tissue
- › reusable or hybrid instrument
- › special durable stainless steel
- › non-breakable jaws
- › double-linked force transfer

OTHER COMMON FEATURES

- › demountable (sheath, insert, handle)
- › 45° handle-top HF-port
- › sheath lengths of 240 mm, 310 mm or 430 mm (no problem with obese patients)
- › 360° rotation, therefore optimal positioning
- › highest durability, even if autoclaved
- › virtually no stock needed, compared to disposable instruments
- › environment-friendly due to its durability
- › compatible with all our inserts for 5 or 10 mm sheaths
- › various handle designs available